




Surface Measurement Systems
World Leader in Sorption Science

Predict Key Cosmetics Material Properties with Advanced Sorption Characterization



Sorption Behavior
& Hygroscopicity



Powder Flow &
Cohesion



Surface Energy &
Wettability



Vapor Pressure of
Fragrances



Absorption &
Release Capacity



Cyclability &
Stability

Our groundbreaking instruments employ Dynamic Vapor Sorption, Inverse Gas Chromatography, & Breakthrough Analysis to provide comprehensive & precise characterization of solid-state materials, empowering research & enabling scientific breakthroughs in pharmaceutical drug development.





Surface Measurement Systems
World Leader in Sorption Science

Empowering Material Research for World's Leading Cosmetics Developers

AVEDA



Johnson & Johnson

L'ORÉAL

TRI Princeton
RESEARCH EDUCATION SERVICE

ESTÉE LAUDER




kenvue


P&G




Our Instruments in Action



Dispersive surface energy distributions of hair samples



Wettability, γ_s^{AB}/γ_s^T profiles for the hair samples.



The superimposed moisture sorption and desorption kinetic results for hair samples; undamaged caucasian virgin (red), bleached

Case Study 616: Determining of Surface Properties of Human Hair by iGC & DVS

The Problem: Hair undergoes surface damage from chemical treatments (e.g., bleaching) and thermal exposure, which alters its physical and chemical behavior. These surface changes affect hydration, adhesion, and product performance, making it essential for the personal care industry to understand and quantify surface energy, wettability, and moisture sorption characteristics of damaged versus untreated hair.

The Techniques:

- Dynamic Vapor Sorption (DVS)
- Inverse Gas Chromatography (iGC)

The Research: Inverse Gas Chromatography (iGC SEA) and Dynamic Vapor Sorption (DVS) were used to measure surface energy heterogeneity, acid-base properties, and water sorption behavior of virgin, thermally treated, and bleached human hair samples.

The Results: Bleached and heat-damaged hair showed increased surface energy, wettability, and water uptake compared to undamaged hair. iGC revealed greater polar interactions in damaged hair due to higher surface heterogeneity. DVS showed significant hysteresis in thermally treated hair, suggesting restricted moisture diffusion. These insights guide formulation of effective hair care treatments.

Scan to Download
Cosmetic Materials
Reference Pack



www.surfacemeasurementsystems.com